Microfluidics based Worm Sorter

Koduri Mani Manohar Prasad¹, John Yan², Andrew A.Bettiol², K.C James Raju¹ 1.University of Hyderabad, Hyderabad, Andhra Pradesh, India 2.National University of Singapore, Singapore

1. Introduction

- The nematode Caenorhabditis Elegans is a widely adopted model organism for studying various neurobiological processes at the molecular and cellular level in vivo.
- With a small, flexible, and continuously moving body, the manipulation of C.elegans becomes a challenging task.
- ■A lab-on-a-chip (LOC) is a device that integrates one or several laboratory functions on a single chip of only millimetres to a few square centimetres in size is employed for screening of C.elegans.



Figure : A. C.elegans

2. Proposed System

■ In particular, Pressure driven Micro Fluid channels are our interest to sort out C.elegans in a channel, based on GFP profiles by integrating two sets of emitting and detection fibers.




Figure: B: Y- Channel

- ■At some pressure flow from Outlet 1 to Inlet or Outlet 2 can be stopped...
- Is there any Numerical Value of Pressure??

4. Simulation Results

■To simulate the behaviour of fluid in a channel Comsol is used




Figure: C: Velocity

Figure: D: Pressure Distribution

in Y Shaped channel




Figure: D: Velocity Vs Pressure

5.Conclusions

- ■A Y-shaped channel is simulated for the optimum numerical pressure variations such that velocity variations in practice is feasible for sorting.
- ■This allows the device to screen C.elegans based on their optical properties without the help of a microscope, an image tracking algorithm and a computer.

3. Theorotical Approach

 By using Symmetry approach, Bernoulli's Equation, and Continuity Equation

we obtained a Control Equation

By solving the control equation

For change in velocity of +10%

variation in pressure is about 37%

For change in velocity of -10%

Variation in Pressure is about 67%

5.References

1. Xavier Casadevall i Solvas et al.

High-throughput age synchronisation of Caenorhabditis elegans Chem. Commun., 2011,47, 9801–9803

2. Kwanghun Chung et al.

Automated on-chip rapid microscopy, phenotyping and sorting of C. elegans

NATURE METHODS | VOL.5 NO.7 | ,637-643, JULY 2008

3.Adela Ben-Yakare et al.

Microfluidics for the analysis of behavior, nerve regeneration, and neural cell biology in C. elegans

Current Opinion in Neurobiology2009,19:561-567

Excerpt from the Proceedings of the 2013 COMSOL Conference in Bangalore





